Microbiology Kirby-Bauer Susceptibility Testing (Modified)

Materials required:

- Misc. equipment and supplies as needed
- Trypticase soy broth (TSB), 1 ml
- Mueller Hinton Agar (MH) 2
- Sterile swab
- Antibiotic disks, as per culture type
- Metric ruler/caliper
- Staphylococcus aureus ATTC 25922 and Escherichia coli ATTC 25923

What is the Clinical Purpose of performing a Susceptibility Test: (1 pt)

Procedure: DAY 1

- 1. Using your inoculating needle touch 3-5 well isolated colonies of the organism to be tested. Use the needle to inoculate a TSB. Incubate at 35° C for one hour. The broth should be noticeably turbid. (A turbidity standard is normally used.)
- 2. Label your plates as directed. Soak a sterile swab in the TSB for a few seconds. Using the spread plate technique, cover the entire MH plate using overlapping streaks. Turn the plate 90 degrees and repeat. Turn the plate 45 degrees and repeat for a third time. The last inoculating step requires running the swab around the inner edge of the plate. The plate must be covered with a uniform lawn of the bacteria to be tested. Repeat this plate inoculation on a second plate.
- Place the antibiotic dispenser on the MH plate. Dispense the antibiotic disks. Place no more than 4 discs on a plate. See reference chart below as to determine which discs to place on each plate. Using the back side of a loop <u>gently</u> tamp the disks onto the media. Good contact is essential for accurate results. If a disk is out of position on the media <u>do not</u> move it. Incubate at 35° C for 18-24 hours.
- 4. **DAY 2:** Observe your MH plate.

a. If there were obvious gaps in the inoculum, what would this indicate; other than the zones of inhibition? $(0.5\ \text{pt})$

b. If there were two or more morphological differences present, what would this indicate? $(0.5\ \text{pt})$

- 5. On the table provided, list the antibiotics and their abbreviations. Using a metric rule or caliper, measure the zone of inhibition in millimeters for each antibiotic. Record your measurements.
- 6. Using the provided chart interpret the results as: S (sensitive), I (intermediate), or R (resistant). Record your results on the table below.

Name: ___

| Susceptibility results: | | S. aureus | | E. coli | |
|-------------------------|-------|-----------|-----------------|---------|------------------------|
| Antibiotic | Abbr. | Zone | Interp (0.5 pt) | Zone | Interp (0.5 pt) |
| Ampicillin | AMP | | | | |
| Chloramphenicol | С | | | | |
| Erythromycin | E | | | | |
| Gentamicin | CN | | | | |
| Penicillin | Р | | | | |
| Streptomycin | S | | | | |
| Tetracycline | TE | | | | |
| Sulfonamides | SF | | | | |

Zone Diameter Interpretive Standards:

| Antimicrobial Agent | Resistant | Intermediate | Susceptible |
|-----------------------|----------------|--------------|----------------|
| Ampicillin (gram neg) | <u><</u> 13 | | <u>></u> 17 |
| Ampicillin (Staph) | <u><</u> 28 | | <u>></u> 29 |
| Chloramphenicol | <u><</u> 12 | 13-17 | <u>></u> 18 |
| Erythromycin | <u><</u> 15 | 14-22 | <u>></u> 23 |
| Gentamicin | <u><</u> 12 | 13-14 | <u>></u> 15 |
| Penicillin (Staph) | <u><</u> 28 | | <u>></u> 29 |
| Streptomycin | <u><</u> 11 | 12-14 | <u>></u> 15 |
| Tetracycline | <u><</u> 14 | 15-18 | <u>></u> 19 |
| Sulfonamides | <u><</u> 12 | | <u>></u> 26 |

Organism Identification (Either Staph or E. coli):

Type your answers to the following questions on a separate sheet of paper and attach.

- 1. Your experimental results should have shown that some of the antibiotics have a broad range of activity against bacteria. Why do these same antibiotics have no activity against eukaryotic cells such as fungi? (1 pt)
- All aspects of the Kirby-Bauer (KB) test are standardized to assure reliability.
 a. What might the consequence be of pouring a plate 2 mm deep instead of 4mm? (1 pt)

b. The plates are supposed to be used within a specific time after their preparation and should be free of visible moisture. What to you surmise as the reason for these directions? (1 pt)

- 3. In clinical applications of the KB test, diluted cultures must be used within thirty minutes. Why is this important? (Hint: look at generation times as explained in chapter 7). (1 pt)
- 4. *E. coli* and *S. aureus* were chosen to represent Gram negative and Gram positive bacteria, respectively. For a given antibiotic, is there a difference in susceptibility between Gram negative and Gram positive bacteria? Explain what the difference(s) is(are) between these two organisms. (1 pt)
- 5. Suppose you do this test on a hypothetical *Staphylococcus* species with the antibiotics penicillin (P-10) and chloramphenicol (C-30). You record zone diameters of 20 mm for the chloramphenicol disc and 25 mm for the penicillin disc. Which antibiotic would be the most

Name: _____

effective against this organism? What does this tell you about comparing zone diameters to each other and the importance of the interpretive chart? (2 pt)